

BIOMÉRIEUX

# AIR IDEAL<sup>®</sup> 3P<sup>®</sup>

I am powerful. I am simple



Your Ally in Advancing Quality

PIONEERING DIAGNOSTICS

*“At bioMérieux we strive to provide the best and most innovative solutions for the biopharmaceutical industries, developed in partnership with experts in Pharmaceutical sciences”*

Skan AG, Johnson Diversey and Health Protection Agency (UK)



# THE SOLUTION FOR AIR SAMPLING ENVIRONMENTS

## Standard of Performances

The AIR IDEAL<sup>®</sup> 3P<sup>®</sup> was third party validated by the Health Protection Agency to meet the stringent requirements of ISO 14698-1 and EN 17141 for the control of clean rooms

### EASE OF USE

- Lightweight (1.2 kg)
- Ergonomic shape and easy to handle
- Audible and visual signals
- Long battery life, short charge time (3 hours)

### SIMPLE INTERFACE

- Programmable
- Rapid set up time, previous airflow program is saved
- Sequential sampling and volume adjustable up to 2000 L

### ROBUSTNESS AND RELIABILITY

- Precise air intake
- Robust with shockproof ABS plastic
- Minimal disruption of unidirectional airflow
- Flow rate and linear velocity: 100L/min and 20m/s
- High autonomy, at least 24 consecutive 1000 liter samples to be collected

### FREE OF CONTAMINATION

- Smooth surface, free of dead angles for ease of cleaning
- Simple, validated air circuit decontamination procedure
  - Disinfectant compatibility
  - Autoclavable sampling grids
  - VHP proof




# 3P® FOR PHARMACEUTICAL PROVEN PERFORMANCES

The AIR IDEAL® 3P® engine, electronic components and seals meet the strict constraints of isolator.

The use of the AIR IDEAL® 3P® in isolators was third party validated by SKAN A.G and LCB and guarantees that the performances are not affected by sequences of isolator decontaminations.

## BIOLOGICAL EFFICIENCY

The bioMérieux AIR IDEAL® 3P® microbial air sampler has been third part validated by the HPA for biological efficiency following the ISO 14698-1 standard.



The air IDEAL® 3P® was shown to be able to collect aerosols of *Staphylococcus epidermidis* with a high efficiency of 92.1%.

This study demonstrates the superior performances of air IDEAL® 3P® on both the strains of reference tested: *B. subtilis* and *S. epidermitis* (indicators of physical and biological efficiency).




## VALIDATION OF THE AIR CIRCUIT DECONTAMINATION EFFICIENCY

The efficiency of the equipment's air circuit decontamination with IPA has been validated following the recommendation of the NF EN 1040 standard.

This procedure allows the reduction of 4.6 log of *P. aeruginosa*.

## PHYSICAL EFFICIENCY

Thirty years of investigation on environmental airborne particles containing bacteria have demonstrated that the range of particle sizes is from 5 to 15µm, with a median size of particles at 13µm.



The result of physical evaluation of AIR IDEAL® 3P® following the ISO 14698-1 standard, demonstrated that the bioMérieux air sampler is highly efficient for collecting particles containing bacterial cells in the range 2.1-14 microns.

Indeed, in this range of particle size, AIR IDEAL® 3P® collects from 85% to 139% of the particles. This study confirmed that air IDEAL® 3P® has a high level of collection of the particles of interest (above 5µm) with an efficiency of collection superior to 100%.



For a fully integrated environmental monitoring solution and a reliable standard of quality, bioMérieux provides a global range of products benefiting from over 50 years of microbiology expertise



## Get the best performance when combined with our 3P® culture media

### AIR IDEAL® 3P®

AIR IDEAL® 3P® equipped with 90 mm diameter sampling grids

**Reference 423154**

### ACCESSORIES

90 mm diameter additional sampling grids

**Reference 96309**

Aluminium tripod with ball joint and adapter

**Reference 96308**

### CULTURE MEDIA

Trypcase soy 3P® irradiated agar with neutralizers and  $\beta$ -lactamase  
**Reference 43287**

Irradiated Trypcase Soy Agar 3P®  
**Reference 43711/43169 and 423725**

Irradiated Trypcase Soy Agar 3P® with Neutralizers  
**Reference 43811/43819 and 423723**

Irradiated Sabouraud Dextrose Agar 3P® with Neutralizers  
**Reference 43814**

