



TEMPO[®]

FOOD QUALITY AUTOMATION

Because you count



PIONEERING DIAGNOSTICS

SIMPLIFY YOUR LIFE

ONLY **1 CARD**


HIGHLY FLEXIBLE!

- MANAGEMENT OF SURGE TESTING NEEDS
- FAST START FOR «LAST MINUTE» SAMPLES
- QUICK TECHNICAL TRAINING


ONLY **3** STEPS

as simple as 1, 2, 3!

- 1** Automated sample preparation of samples



- 2** Incubation in your standard incubators




- 3** Automated reading and recording



«WITH **TEMPO®**,
I CAN **TEST ON DEMAND**»

REDUCE

- WASTE
Protect the environment and be lean
- INCUBATOR SPACE
- SUBJECTIVITY



RESULTS YOU CAN COUNT ON

GLOBALLY RECOGNIZED

- AFNOR & AOAC VALIDATED
- NO CROSS CONTAMINATION: SECURED PROTOCOL
- LONG SHELF LIFE REAGENTS

DATA INTEGRITY

KEEP YOUR DATA WITHIN COMPLIANCE STANDARDS

- FULL TRACEABILITY - INTEGRATED CONTROL PLAN & BARCODE
- SECURED DATA - FROM SAMPLE ID TO AUTOMATED RESULTS
- NO TRANSCRIPTION ERRORS

— bioMérieux —
CONNECT-UP™

MINIMIZE ADMINISTRATIVE TASKS,
GAIN ON VALUE-ADDED ACTIONS!

TEMPO® parameters


TEMPO parameters	NEXT DAY SOLUTION	AOAC	ISO
411113 AC Aerobic mesophilic total flora	✓	✓	✓
80002 STA Staphylococcus (coag+)	✓	✓	✓
80003 EB Enterobacteriaceae	✓	✓	✓
80004 EC Escherishia coli	✓	✓	✓
80006 TC Coliforms	✓		✓
80044 CC Coliforms	✓	✓	
80106 BC Bacillus cereus	✓	✓	✓
80001 YM Yeast & Mold	72-76 hr	✓	
80071 LAB Lactic acid bacteria	40-48 hr		
421509 CAM Campylobacter	48 hr	✓	

GET THE RIGHT SOLUTION

SAMPLE PREPARATION



PRODUCTIVITY



SAFETY



IDENTIFICATION

