



TEMPO®

List of validated methods

Reference	Target	Method	ISO 16140-2 (2016)			AOAC RI	AOAC OMA	Health Canada	Other validations
			Reference method	Certificate #	Expiration	Certificate #	Certificate #	Certificate #	Certificate #
80106	Enumeration <i>B. cereus</i>	TEMPO® BC	ISO 7932	MICROVAL 2014LR47	2027 12	PTM 071401			
80004	β -glucuronidase-positive <i>Escherichia coli</i>	TEMPO® EC	ISO 16649-2	AFNOR BIO 12/13-02/25	2025 02	PTM 080603	OMA 2009.02	MFLP-107 (2018)	India (FSSAI), Rapid Analytical Food Testing (RAFT): Direct approval (exp. 10.2025) Vietnam Standards (TCVN): Official approval Australia (Department of Agriculture, Water and the Environment): Official approval AOAC 2009.02
411113	Colony-count at 30 °C	TEMPO® AC	ISO 4833-1	AFNOR BIO 12/35-05/13	2025 05	PTM 121204		MFLP-105 (2018)	India (FSSAI), Rapid Analytical Food Testing (RAFT): Direct approval Vietnam Standards (TCVN): Official approval Australia (Department of Agriculture, Water and the Environment): Official approval AOAC 2008.10
80006	Coliforms (37°C)	TEMPO® TC	ISO 4832	AFNOR BIO 12/17-12/05	2025 12				
80002	Coagulase-positive Staphylococci	TEMPO® STA	ISO 6888-2	AFNOR BIO 12/28-04/10	2026 04	PTM 120901		MFLP-109 (2020)	India (FSSAI), Rapid Analytical Food Testing (RAFT): Direct approval (exp. 02.2027)
80003	Enterobacteriaceae	TEMPO® EB	ISO 21528-2	AFNOR BIO 12/21-12/06	2026 12	PTM 050801			India (FSSAI), Rapid Analytical Food Testing (RAFT): Direct approval (exp. 02.2027)
80044	Coliforms	TEMPO® CC				PTM 060702		MFLP-106 (2018)	
80001	Yeasts and molds	TEMPO® YM				PTM 041001		MFLP-108 (2018)	
421509	<i>Campylobacter</i> spp.	TEMPO® CAM	ISO 10272-2	AFNOR BIO 12/43-04/20	2028 04	PTM 112002			

bioMérieux S.A. • 69280 Marcy l'Étoile • France • Tel: +33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90
www.biomerieux.com